

DISCOVER

MK-dent®
GERMANY



HIGH- & LOW SPEED HANDPIECES | PROPHYLAXIS | SCALER
MOTORS | QUICK CONNECTORS | TUBINGS | SPARE PARTS | TOOLS

www.MK-dent.com



THE PRIME LINE SERIES

MAX. 55dB(A): QUIET LIKE NO OTHER!

The Prime Line spindle has a two-part bur guide with a carbide shell for the highest level of concentricity and resistance to bending. With over 40 N of retention force, the spindle guarantees that the bur is securely held. The vibration-free running results in a sound level of less than 55 dB. This makes the Prime Line quiet like no other (HP21KL).

FEATURES

- Prime Grip
- Premium spindle with carbide-insert
- Pressure reduction system (only handpieces with KaVo Multiflex® connection)
- Titanium coating
- Ceramic bearings
- Pushbutton system
- Glass rod light conductor
- Non-return valves

HIGH SPEED HANDPIECES



595

HP21KL
High Speed Handpiece with light
Standard head
25W
4-hole spray
KaVo Multiflex® connection

595

HP21NL
High Speed Handpiece with light
Standard head
25W
4-hole spray
NSK Phatelus® connection



595

HP22KL
High Speed Handpiece with light
Small head
18W
3-hole spray
KaVo Multiflex® connection

595

HP22NL
High Speed Handpiece with light
Small head
18W
3-hole spray
NSK Phatelus® connection



THE CLASSIC LINE SERIES

FEEL THE ADVANTAGE!

One result of our research is the newly designed grip for the Classic Line. The Classic Grip fits perfectly in your hand. Tangibly improved ergonomics and higher stability are the result. A high resistance to wearout and corrosion is ensured through the high quality coating of the surface.

FEATURES

- Classic Grip
- Premium spindle with carbide-insert
- Chrome coating
- Ceramic bearings
- Pushbutton system
- Glass rod light conductor

HIGH SPEED HANDPIECES WITH LIGHT



495

HC20KL - Power Edition
High Speed Handpiece with light
Standard head
28W
4-hole spray
KaVo Multiflex® connection

28W!



495

HC22KL
High Speed Handpiece with light
Small head
18W
3-hole spray
KaVo Multiflex® connection



395

HC20K - Power Edition
High Speed Handpiece
Standard head
28W
4-hole spray
KaVo Multiflex® connection

28W!



395

HC22K
High Speed Handpiece
Small head
18W
3-hole spray
KaVo Multiflex® connection

HC21K
High Speed Handpiece
Standard head
23W
4-hole spray
KaVo Multiflex® connection

395

THE ECO LINE SERIES

HIGH FLEXIBILITY, LOW PRICE!

FEATURES

- Smooth Grip
- Ceramic bearings
- Spindle with steel-insert
- Pushbutton system
- Chrome coating
- Glass rod light conductor



HIGH SPEED HANDPIECES WITH MULTIFLEX® CONNECTION



HE21KL **395**
High Speed Handpiece with light
Standard head
21W
4-hole spray
KaVo Multiflex® connection

HE21NL **395**
High Speed Handpiece with light
Standard head
21W
4-hole spray
NSK Phatelus/Mach® connection



HE22KL **395**
High Speed Handpiece with light
Small head
16W
3-hole spray
KaVo Multiflex® connection

HE22NL **395**
High Speed Handpiece with light
Small head
16W
3-hole spray
NSK Phatelus/Mach® connection



HE21K **295**
High Speed Handpiece
Standard head
21W
4-hole spray
KaVo Multiflex® connection

HE21N **295**
High Speed Handpiece
Standard head
21W
4-hole spray
NSK Phatelus/Mach® connection



HE22K **295**
High Speed Handpiece
Small head
16W
3-hole spray
KaVo Multiflex® connection

HE22N **295**
High Speed Handpiece
Small head
16W
3-hole spray
NSK Phatelus/Mach® connection

HIGH SPEED HANDPIECES WITH TUBING CONNECTION



HE11 249
High Speed Handpiece
Standard head
18W
1-hole spray
2-hole connection



HE12 249
High Speed Handpiece
Standard head
18W
1-hole spray
4-hole connection



HE14 249
High Speed Handpiece
Mini head
15W
1-hole spray
2-hole connection



**NIET VERKRIJGBAAR
IN BELGIË**
HE15L
High Speed Handpiece with light
Standard head
18W
4-hole spray
5-hole connection



**NIET VERKRIJGBAAR
IN BELGIË**
HE16L
High Speed Handpiece with light
Small head
16W
3-hole spray
5-hole connection



HE17 295
High Speed Handpiece
Standard head
18W
4-hole spray
4-hole connection



HE18 295
High Speed Handpiece
Small head
16W
3-hole spray
4-hole connection



THE BASIC LINE SERIES

PROVEN TECHNOLOGY, SAFE GRIP!

FEATURES

- Basic Grip
- Ceramic bearings
- Spindle with steel-insert
- Pushbutton system
- Chrome coating
- Glass rod light conductor



HIGH SPEED HANDPIECES WITH MULTIFLEX® CONNECTION



HB21KL **395**
High Speed Handpiece with light
Standard head
21W
4-hole spray
KaVo Multiflex® connection



HB23KL **395**
High Speed Handpiece with light
Mini head
15W
1-hole spray
KaVo Multiflex® connection



HB21K **295**
High Speed Handpiece
Standard head
21W
4-hole spray
KaVo Multiflex® connection



HB23K **295**
High Speed Handpiece
Mini head
15W
1-hole spray
KaVo Multiflex® connection

HIGH SPEED HANDPIECES WITH TUBING CONNECTION



HB12 **249**
High Speed Handpiece
Standard head
18W
1-hole spray
4-hole connection



HB14 **249**
High Speed Handpiece
Mini head
15W
1-hole spray
4-hole connection





THE PRIME LS LINE SERIES

PRIME GRIP: OPERATING CONVENIENCE!

The combination of a brushed stainless steel housing and a titanium coating offers a perfect grip, protection against corrosion and maximum scratch resistance.

Due to the ideal angle combination, you can reach hard to reach preparation areas for maximum comfort during treatment.

CONTRA ANGLE HEADS

FEATURES

- KaVo® compatible
- Titanium coating
- Push button



198
LPH11
Contra Angle Head
1:1, blue
Max. 40.000 rpm
For CA burs (Ø2.35 mm)



239
LPH16
Contra Angle Head
1:1, blue
Max. 40.000 rpm
For FG burs (Ø1.6 mm)



198
LPH21
Contra Angle Head
2:1, green
Max. 20.000 rpm
For CA burs (Ø2.35 mm)

STRAIGHT- & CONTRA ANGLE HANDPIECES

FEATURES

- Prime Grip
- With light
- Intramatic® connection
- Stainless steel case with titanium coating



450
LP01L
Straight Handpiece with light
1:1, blue
With internal water
Max. 40.000 rpm
For HP burs (Ø2.35 mm)



350
LP11L
Contra Angle Body with light
1:1, blue
With internal water
Max. 40.000 rpm
For MK-dent and KaVo® heads



450
LP27L
Contra Angle Body with light
2,7:1, green
With internal water
Max. 14.814 rpm
For MK-dent and KaVo® heads



495
LP74L
Contra Angle Body with light
7,4:1, green
With internal water
Max. 5.406 rpm
For MK-dent and KaVo® heads



THE CLASSIC LS LINE SERIES

CLASSIC GRIP: FEEL THE ADVANTAGE!

The Classic Grip fits perfectly in your hand. Tangibly improved ergonomics and higher stability are the result. A high resistance to wear and corrosion is ensured through the high quality chrome coating of the surface and the brass housing.

CONTRA ANGLE HEADS

FEATURES

- KaVo® compatible
- Chrome coating
- Push button



LCH11 **185**
Contra Angle Head
1:1, blue
Max. 40.000 rpm
For CA burs (Ø 2.35 mm)



LCH16 **225**
Contra Angle Head
1:1, blue
Max. 40.000 rpm
For FG burs (Ø 1.6 mm)



LCH21 **185**
Contra Angle Head
2:1, green
Max. 20.000 rpm
For CA burs (Ø 2.35 mm)

STRAIGHT- & CONTRA ANGLE HANDPIECES

FEATURES

- Classic Grip
- Available with or without light
- Intramatic® connection
- Chrome coating



LC01L / LC01 **395 / 295**
Straight Handpiece with / without light
1:1, blue
With internal water
Max. 40.000 rpm
For HP burs (Ø 2.35 mm)



LC11L / LC11 **295 / 250**
Contra Angle Body with / without light
1:1, blue
With internal water
Max. 40.000 rpm
For MK-dent and KaVo® heads



LC27L / LC27 **395 / 295**
Contra Angle Body with / without light
2,7:1, green
With internal water
Max. 14.814 rpm
For MK-dent and KaVo® heads



LC74L / LC74 **395 / 345**
Contra Angle Body with / without light
7,4:1, green
With internal water
Max. 5.406 rpm
For MK-dent and KaVo® heads



THE ECO LS LINE SERIES SOLID PERFORMANCE!

The Eco Line Contra Angle Handpieces with built-in contra angle head offer high-tech precision and convincing reliability for a nice price.

The wide spectrum of transmissions provides the right handpiece for each purpose.

STRAIGHT- & CONTRA ANGLE HANDPIECES

FEATURES

- Smooth Grip
- Available with or without light
- Intramatic® connection
- Chrome coating



LE01L / LE01 **350 / 250**
Straight Handpiece with / without light
 1:1, blue
 With internal water
 Max. 40.000 rpm
 For HP burs (Ø 2.35 mm)



LE11L / LE11 **395 / 295**
Contra Angle Handpiece with / without light
 1:1, blue
 With internal water
 Max. 40.000 rpm
 For CA burs (Ø 2.35 mm)



LE16L / LE16 **425 / 325**
Contra Angle Handpiece with / without light
 1:1, blue
 With internal water
 Max. 40.000 rpm
 For FG burs (Ø 1.6 mm)



LE21L / LE21 **395 / 295**
Contra Angle Handpiece with / without light
 2:1, green
 With internal water
 Max. 20.000 rpm
 For CA burs (Ø 2.35 mm)



LE27L / LE27 **450 / 350**
Contra Angle Handpiece with / without light
 2,7:1, green
 With internal water
 Max. 14.814 rpm
 For CA burs (Ø 2.35 mm)



LE54L / LE54 **450 / 350**
Contra Angle Handpiece with / without light
 5,4:1, green
 With internal water
 Max. 7.407 rpm
 For CA burs (Ø 2.35 mm)



LE74L / LE74 **495 / 395**
Contra Angle Handpiece with / without light
 7,4:1, green
 With internal water
 Max. 5.406 rpm
 For CA burs (Ø 2.35 mm)



LE148L / LE148 **495 / 395**
Contra Angle Handpiece with / without light
 14,8:1, green
 With internal water
 Max. 2.703 rpm
 For CA burs (Ø 2.35 mm)



THE BASIC LS LINE SERIES

PROVEN TECHNOLOGY, SAFE GRIP!

The Basic Grip guarantees a safe and comfortable feeling at work – in all situations. The chrome coating is particularly resistant and easy to clean.

Also, on a selection of Basic Line straight- and contra angle handpieces, you have the choice between external water or no water supply.

STRAIGHT HANDPIECES

FEATURES

- Basic Grip
- Intramatic® connection
- Chrome coating



LB01 **195**
Straight Handpiece
 1:1, blue
 Max. 20.000 rpm
 For HP burs (Ø 2.35 mm)
 Intramatic® connection



LB02 **195**
Straight Handpiece
 1:1, blue
 External water via clip-on hose
 Max. 20.000 rpm
 For HP burs (Ø 2.35 mm) and disposable prophylaxis heads
 Intramatic® connection

CONTRA ANGLE HEADS



LBH11 95
Contra Angle Head
1:1, blue
Max. 20.000 rpm
Pushbutton system
For CA burs (Ø 2.35 mm)



LBH11P 59
Contra Angle Head
1:1, blue
Max. 5.000 rpm
For snap-on prophylaxis cups



LBH11C 59
Contra Angle Head
1:1, blue
Max. 20.000 rpm
Latch type system
For CA burs (Ø 2.35 mm)

CONTRA ANGLE HANDPIECES



LB11 85
Contra Angle Bottom Part
1:1, blue
Max. 20.000 rpm
For MK-dent LBH11/11P/11C and
NSK® EX series heads
Intramatic® connection



LB41 139
Contra Angle Bottom Part
4:1, green
Max. 5.000 rpm
For MK-dent LBH11/11P/11C and
NSK® EX series heads
Intramatic® connection



LB11/11 195
Contra Angle Handpiece
1:1, blue
External water via clip-on hose
Max. 20.000 rpm
For CA burs (Ø 2.35 mm)
Intramatic® connection



LB11/11C 149
Contra Angle Handpiece
1:1, blue
External water via clip-on hose
Max. 20.000 rpm
For CA burs (Ø 2.35 mm)
Intramatic® connection



LB41/11 195
Contra Angle Handpiece
4:1, green
External water via clip-on hose
Max. 5.000 rpm
For CA burs (Ø 2.35 mm)
Intramatic® connection



LB41/11P 195
Contra Angle Handpiece
4:1, green
Max. 5.000 rpm
For snap-on prophylaxis cups
Intramatic® connection

THE PROPHY LINE SERIES

LIGHT, EFFICIENT, UNIVERSAL!






FEATURES

- Ergonomically shaped
- Light weight and compact
- 360° rotatable outer handle RotoGrip®-system for simple and fatigue proof operation



Air Flow Handpiece

Available for the following connections:

-  KaVo® (PR1011K)
-  W&H® (PR1011W)
-  Sirona® (PR1011S)
-  NSK® (PR1011N)
-  Bien-Air® (PR1011B)

595



Prophylaxis Cleaning Powder

For Propphy Line and other air flow handpieces

Extra fine grain size: 12 - 14 µm



Available in the following flavors:

- Mint flavor (PP1011)
- Lemon flavor (PP1012)
- Strawberry flavor (PP1013)

18,50

ULTRASONIC SCALER

FEATURES

- Ergonomic shape for excellent grip
- Lean steel sleeve for a long durability
- Painfree treatments for the patient through linear movements of the tip
- 6-LED light ring for optimum lighting
- Available for the following connections:  EMS®  Satelec®



SC21E **195**

Scaler Handpiece

24 - 33 kHz

Steel funnel cover

EMS® connection, type EN-041



SC21EL **295**

Scaler Handpiece with light

24 - 33 kHz

Steel funnel cover

EMS® connection, type EN-060A



SC21S **195**

Scaler Handpiece

24 - 33 kHz

Steel funnel cover

Acteon/Satelec® connection,
type F 12 281



SC21SL **295**

Scaler Handpiece with light

24 - 33 kHz

Steel funnel cover

Acteon/Satelec® connection,
type F 12 609

THE INNOVATIVE DLC COATING TECHNIQUE: DIAMOND-LIKE CARBON!



FEATURES

- High quality DLC coating
- Extreme hard and wear resistant, for an up to 40% longer lifetime
- Optimum view through a reflection-poor surface
- Anti-allergenic



Art. No. for EMS® / Satelec®	ST01E / ST01S	ST01ED / ST01SD	ST02E / ST02S	ST02ED / ST02SD
Application area	Scaling	Scaling	Scaling	Scaling
Coating/Material	Stainless steel	DLC	Stainless steel	DLC



Art. No. for EMS® / Satelec®	ST03E / ST03S	ST03ED / ST03SD	ST04E / ST04S	ST04ED / ST04SD
Application area	Scaling	Scaling	Scaling	Scaling
Coating/Material	Stainless steel	DLC	Stainless steel	DLC



Art. No. for EMS® / Satelec®	ST05E / ST05S	ST06E / ST06S	ST06ED / ST06SD	ST07E / ST07S
Application area	Endodontics	Periodontal	Periodontal	Periodontal
Coating/Material	Stainless steel	Stainless steel	DLC	Stainless steel



49



74



49



49

Art. No. for EMS® / Satelec®	ST08E / ST08S	ST08ED / ST08SD	ST09E / ST09S	ST10E / ST10S
Application area	Periodontal	Periodontal	Periodontal	Periodontal
Coating/Material	Stainless steel	DLC	Stainless steel	Stainless steel




74



49



45



45

Art. No. for EMS® / Satelec®	ST10ED / ST10SD	ST11E / ST11S	ST12E / ST12S	ST13E / ST13S
Application area	Periodontal	Periodontal	Endodontics	Endodontics
Coating/Material	DLC	Stainless steel	Stainless steel	Stainless steel



45



49



49



49

Art. No. for EMS® / Satelec®	ST14E / ST14S	ST15E / ST15S	ST16E / ST16S	ST17E / ST17S
Application area	Endodontics	Endodontics	Endodontics	Endodontics
Coating/Material	Stainless steel	Stainless steel	Stainless steel	Stainless steel



49



54



54



54

Art. No. for EMS® / Satelec®	ST18E / ST18S	ST19E / ST19S	ST20E / ST20S	ST21E / ST21S
Application area	Endodontics	Cavity preparation	Cavity preparation	Cavity preparation
Coating/Material	Stainless steel	Stainless steel	Stainless steel	Stainless steel



54



54

Art. No. for EMS® / Satelec®	ST22E / ST22S	ST23E / ST23S
Application area	Cavity preparation	Cavity preparation
Coating/Material	Stainless steel	Stainless steel



19

SC21KEY Torque Key for insertion and storage of scaler tips

- Combination of torque key and key holder for neat storage and sterilization of scaler tips
- No contact required during handling of the instrument - therefore no risk of injury or contamination
- A defined torque of 0,3 - 0,4 Nm prevents an overwinding of the thread
- Autoclavable at up to 135 °C

QUICK CONNECTORS FOR HIGH SPEED HANDPIECES



6-hole Quick Connector with spray regulation system

MK-dent & KaVo Multiflex® connection

Available in the following versions:

QC6016KW with LED 198

QC5016KW with Xenon 179



6-hole Quick Connector

MK-dent & KaVo Multiflex® connection

Available in the following versions:

QC6016K with LED 169

QC5016K with Xenon 149



Quick Connector

MK-dent & KaVo Multiflex® connection

Available in the following versions:

QC4014K without light (4-hole) 89

QC4012K without light (2-hole) 95



6-hole Quick Connector with spray regulation system

Sirona® connection

Available in the following versions:

QC6016SW with LED 198

QC5016SW with Xenon 179



Quick Connector

W&H Roto Quick® connection

Available in the following versions:

QC5016W with power supply (6-hole) 149

QC4014W without power supply (4-hole) 89



Quick Connector

NSK MachLite® connection

Available in the following versions:

QC6016NT with LED (6-hole) 129

QC5016NT with Xenon (6-hole) 98

QC4014NT without light (4-hole) 89

FOR ALL STRAIGHT & CONTRA ANGLE HANDPIECES: POWERFUL DRIVE!



FEATURES

- Electric Motor **EM1190K** with LED
- Max. 40.000 rpm
- Spray regulation system
- Guided interior separate coolants
- Short Intramatic® 2 connection, suitable for all common instruments
- Compatible with KaVo® 190, 196 and K200 connections

695

AIRMOTORS WITH INTERNAL WATER SUPPLY



FEATURES

- 4-hole Airmotor **AM1004**
- Max. 20.000 rpm
- Internal water supply
- Infinitely adjustable right and left drive
- Guided interior separate coolants
- Standard Intramatic® connection, suitable for all common instruments

295

AIRMOTORS WITH EXTERNAL WATER SUPPLY



AM0012 **195**
2-hole Airmotor with external water supply
Max. 20.000 rpm
Standard Intramatic® connection



AM0014 **195**
4-hole Airmotor with external water supply
Max. 20.000 rpm
Standard Intramatic® connection

AIRMOTORS WITHOUT WATER SUPPLY



AM0002 **195**
2-hole Airmotor without water supply
Max. 20.000 rpm
Standard Intramatic® connection



AM0004 **195**
4-hole Airmotor without water supply
Max. 20.000 rpm
Standard Intramatic® connection

REPLACEMENT CARTRIDGES: BREATHAKING PERFORMANCE!



CHOOSE BETWEEN THREE QUALITIES:

- Premium Cartridge - Carbide insert with ceramic bearings
- Standard Cartridge - Steel insert with ceramic bearings
- Basic Cartridge - Steel insert with steel bearings

Your Turbine	Suitable Premium Cartridge
MK-dent	
HB12	TC2021 189
HB14	TC2033 189
HB15L	TC2021 189
HB16L	TC2032 189
HB21K/KL	TC2034 189
HB23K/KL	TC2033 189
HC20K/KL	TC2038 189
HC21K/KL	TC2038 189
HC22K/KL	TC2039 189
HE11	TC2021 189
HE12	TC2021 189
HE14	TC2035 189
HE15L	TC2021 189
HE16L	TC2035 189
HE17	TC2034 189
HE18	TC2035 189
HE21K/KL/N/NL	TC2034 189
HE22K/KL/N/NL	TC2035 189
HP21KL/NL	TC2041 189
HP22KL/NL	TC2042 189
KaVo®	
4500B/BR	TC2500BK 189
5000/B/BY	TC2500BK 189
6000B	TC2600BK 189
630B/C	TC2640BK 189
632D/E/EA	TC2642BK 189
635B	TC2635BK 189
635BR	TC2637BK 189
636CP/P	TC2636BK 189
637B/C	TC2637BK 189
639B/C	TC2642BK 189
640B/C	TC2640BK 189
642B/C	TC2642BK 189
645B/C	TC2642BK 189
646B/C	TC2642BK 189

Your Turbine	Suitable Premium Cartridge	
KaVo®		
647B	TC2647BK	189
649B	TC2647BK	189
6500B	TC2650US	189
6500BR	TC2800BK	189
650B	TC2650BK	189
655B/C	TC2655BK	189
659BR	TC2655BK	189
660B/C	TC2655BK	189
7000B/C	TC2700BK	189
8000B/C/BM/BS/BN	TC2800BK	189
8700	TC2870BK	189
9000	TC2900BK	189
E675L	TC2677BK	189
E677L	TC2677BK	189
E679L	TC2680BK	189
E680L	TC2680BK	189
W&H®		
TA96	TC2012W	189
TA97 C LED	TC2097WT	189
TA98	TC2011W	189
TA98 C LED	TC2098WT	189
TK100	TC2100WK	189
TK98	TC2098WK	189
Sirona®		
T1 Control TC3	TC2011S1	189
T2/T3 Racer	TC2011S2R	189
NSK®		
X600	TC2600NX	189
X600L	TC2600NX	189
X700	TC2700NX	189
X700L	TC2700NX	189
Bien-Air®		
Black Pearl	TC2011BBP	189
Bora	TC2011BB	189
Boralina	TC2011BB	189
Prestige	TC2011BP	189
B.A. International®		
Ultimate 675	TC2011S2R	189
Ultimate 685	TC2011S2R	189
Fona®		
6080MR	TC2011S2R	189
8060SR	TC2011S2R	189

LED: ULTRA HIGH BRILLIANCE!



FEATURES

- Excellent brilliance, no glare or haloes
- 5500 K daylight color temperature
- Clear color contrast, therefore very eye-friendly
- Flexible, user-friendly and economic
- 15 times longer durability compared to other halogen and XENON illuminants
- Available for the following brands:

K KaVo® W W&H® S Sirona® NT NSK® B Bien-Air® CA Castellini® ST Star® F Faro® Y Yoshida® M Midwest®

ALLEMAAL 39

BU8012	For MK-dent and KaVo® quick connectors
BU8012KM	For MK-dent AM1016 and KaVo® motors
BU8012EM	For MK-dent EM1190K
BU8012SB	For Sirona® quick connectors, SPRAYVIT L tubings and Sirosonic L scaler
BU8012SG	For Sirona® motors
BU8012BA	For Bien-Air® motors (not for rotating current)
BU8012BA/V2*	For Bien-Air® motors (not for rotating current)
BU8012BAK	For Bien-Air® quick connectors and high speed handpieces, Midwest® 6-Pin quick connectors
BU8012WH	For W&H® high speed handpieces
BU8012WH/V2*	For W&H® high speed handpieces
BU8012WLT	For W&H® airmotors LT25
BU8012WEA	For W&H® motors 40LT, 50LT and A51LT
BU8012WRA	For W&H® quick connectors RA24
BU8012WRA/V2*	For W&H® quick connectors RA24
BU8012N	For NSK® quick connectors
BU8012NT	For NSK® motors TIM40 und TI205L
BU8012NNL	For NSK® NL motors
BU8012NM	For NSK® motors M40 LED, M205LG, NBS, X205L and quick connectors: KCL, SCL, PTL LED
BU8012C	For Castellini® high speed handpieces
BU8012C/V2*	For Castellini® high speed handpieces
BU8012ST	For Star® quick connectors
BU8012ST/V2*	For Star® quick connectors
BU8012F	For Faro® high speed handpieces
BU8012Y	For Yoshida® quick connectors

* V2 - LEDs with reversed polarity

XENON:

HIGH LUMINOSITY, FAVOURABLE PRICE!



FEATURES

- High luminosity
- Precise illuminance
- Long durability

• Available for the following brands: **K** KaVo® **W** W&H® **S** Sirona® **NT** NSK® **B** Bien-Air® **ST** Star®

ALLEMAAL 12,50

BU7012	For MK-dent and KaVo® quick connectors and motors
BU7012SB	For Sirona® quick connectors, SPRAYVIT L tubings and Sirosonic L scaler
BU7012SG	For Sirona® motors
BU7012WH	For W&H® high speed handpieces
BU7012BA	For Bien-Air® motors
BU7012BAK	For Bien-Air® high speed handpieces
BU7012N	For NSK® quick connectors
BU7012ST	For Star® quick connectors



PREMIUM SERVICE OIL: OPTIMUM CARE!



Full synthetic Premium Service Oil, 500ml

For maintenance of high and low speed handpieces, airmotors and air flow handpieces. With mint scent.

Available in the following versions:

LU1011 suitable for a large selection of lubrication tools for many brands

LU1022 suitable for the KaVo® Quattrocare 2 system

LU1031 suitable for the W&H® Assistina and other maintenance units

18,90
18,90
35

ALUMINIUM LUBRICATION TOOLS

	4,5 	4,5 	4,5 	4,5 
	LT1011	LT1012	LT1013	LT1014
Suitable for	High speed handpieces with MK-dent connection	High speed handpieces with KaVo® connection	Low speed handpieces with Intramatic® connection	FG and PB spindles, Ø 1.6 mm
	6,5 	4,5 	4,5 	4,5 
	LT1015	LT1017	LT1018	LT1019
Suitable for	4-hole and 5-hole connections	High speed handpieces with Sirona® connection	Handpieces with W&H® connection	Handpieces with NSK® connection (CHQD)
	4,5 	4,5 	4,5 	
	LT1020	LT1021	LT1022	
Suitable for	Handpieces with NSK® connection (Ti-Max)	High speed handpieces with Star® connection	Handpieces with Yoshida® connection	

KEYS

22



KE1011KT Back Cap Key

Suitable for KaVo® Super Torque/Bella Torque high speed handpieces with 2-hole back caps. With attachment for torque wrench RT2025 (1,0 Nm).

22



KE1022KT Spray Cover Key

Suitable for MK-dent & KaVo® high speed handpieces with standard head. With attachment for torque wrench RT2025 (0,3 Nm).

22



KE1023T Back Cap Key

Suitable for MK-dent high speed handpieces with small head. With attachment for torque wrench RT2025 (1,0 Nm).

22



KE1637T Back Cap Key

Suitable for KaVo® Mira 635/Bella Torque 637 high speed handpieces. With attachment for torque wrench RT2025 (1,0 Nm).

22



KE1650T Back Cap Key

Suitable for MK-dent high speed handpieces with standard head and 3-/4-hole spray and KaVo® Super Torque/Gentle Force high speed handpieces. With attachment for torque wrench RT2025 (1,0 Nm).

14,50



KE1922T Spray Cover Key

Suitable for MK-dent Prime Line high speed handpieces with standard head. With attachment for torque wrench RT2025 (0,35 Nm).

22



KE2011BT Key

Suitable for low speed handpieces (SP2288 & SP1021). With attachment for torque wrench RT2026 (1,8 Nm).

22



KE2011T Back Cap Key

Suitable for low speed handpieces. With attachment for torque wrench RT2025 (1,0 Nm)

MK-dent[®]
GERMANY



MK-dent[®]
GERMANY

Eichenweg 7b
D-22941 Bargteheide, Germany

Phone: +49 (0)4532 / 400 49-0

Fax: +49 (0)4532 / 400 49-10

Email: info@mk-dent.com

Your distributor

www.MK-dent.com